

Manual: Aunis - Nanonis Scripting Interface

Version: 0.41 (28.08.2025)

Taner Esat



Contents

Manual: Aunis - Nanonis Scripting Interface	1
1) Interface Overview	2
2) Scripting Commands	4
General.....	4
Feedback.....	4
Bias and Current.....	4
Scan	5
Spectroscopy and Lock-In	5
Atom tracking.....	6
Tip position	7
Drift compensation	8
3) Tip Manipulation	9
4) External TCP Interface	10
5) Adding New Functions to the Command System.....	11
Structure of a Registry Entry	11
Examples	12
Adding a Remote Command	13

1) Interface Overview

1. Status

Displays whether the connection to the Nanonis software is established.

Details:

- If connected, the current **bias voltage** and **setpoint current** are shown.
- **Feedback control** can be toggled **ON/OFF** by pressing the button next to the *Feedback* label.
- To establish a connection, press the **Connect** button.
- TCP parameters can be configured under **Settings**.

2. Scripting

Provides a text area for entering and executing script commands.

Usage:

- Enter one or more commands (see [Scripting Commands](#) section).
- The system performs a syntax check before execution to ensure that all commands are valid.
- Press **Run** to execute commands sequentially.
- Press **Stop** to halt execution (note: the current command runs to completion).
- Use the **File** menu to save or load scripts.
- Use the **Autocomplete** feature:
 - Suggested commands appear while typing.
 - Use **arrow keys** and **Tab** to select a suggestion.

3. Tip Manipulation

Allows manual positioning of the tip along the **x**, **y**, and **z** directions.


- **Reference:** For detailed instructions, see section [Tip Manipulation](#).

4. External

Displays an overview of the external TCP interfaces.

5. Settings

Provides configuration options for the TCP connection with the Nanonis software.



2) Scripting Commands

General

Command	Parameter: (Unit)	Description and Example
loop [number] ... end	[number]: (none)	Executes the commands between “loop” and “end” [number] times. Nested loops are also possible. Example: loop 10 ... end
wait [time]	[time]: (s)	Waits [time] seconds before executing the next command. Example: wait 10

Feedback

Command	Parameter: (Unit)	Description and Example
fb.Get	-	Returns the status of the Z-Controller. 1: Feedback loop ON 0: Feedback loop OFF Example: fb.Get
fb.Set [status]	[status]: (none)	Sets the status of the Z-Controller to [status]. [status]: 1 – Feedback loop ON [status]: 0 – Feedback loop OFF Example: fb.Set 1

Bias and Current

Command	Parameter: (Unit)	Description and Example
bias.Get	-	Returns the bias voltage in Volts. Example: bias.Get
bias.Set [bias voltage]	[bias voltage]: (V)	Sets the bias voltage to [bias voltage]. Example: bias.Set 0.5

bias.Add [bias voltage]	[bias voltage]: (V)	Adds [bias voltage] to the current bias voltage. Example: bias.Add 0.1
current.Get	-	Returns the setpoint current in Ampere. Example: current.Get 0.1
current.Set [current]	[current]: (A)	Sets the setpoint current to [current] Example: current.Set 200e-12
current.Add [current]	[current]: (A)	Adds [current] to the setpoint current. Example: current.Add 100e-12

Scan

Command	Parameter: (Unit)	Description and Example
scan.Start	-	Starts a new scan of an image. Note: Next command in the list is executed directly, it does not wait until the scan is finished. Example: scan.Start
scan.Wait	-	Waits until the current scan is finished. Next command in the list is not executed until the scan is finished. Example: scan.Wait

Spectroscopy and Lock-In

Command	Parameter: (Unit)	Description and Example
biasSpec.Start	-	Performs bias spectroscopy. Example: biasSpec.Start
biasSpec.LimitsGet	-	Returns the start and end values for the bias spectroscopy. Example: biasSpec.LimitsGet

biasSpec.LimitsSet [start] [end]	[start]: (V) [end]: (V)	Sets the start and end values for bias spectroscopy to [start] and [end] respectively. Example: biasSpec.LimitsSet -250e-3 50e-3
lockin.FreqGet	-	Returns the frequency of the Lock-In (demod. number 1). Example: lockin.FreqGet
lockin.FreqSet [frequency]	[frequency]: (Hz)	Sets the frequency of the Lock-In (demod. number 1) to [frequency]. Example: lockin.FreqSet 187.0
lockin.AmplGet	-	Returns the amplitude of the Lock-In (demod. number 1). Example: lockin.AmplGet
lockin.AmplSet [amplitude]	[amplitude]: (V)	Sets the amplitude of the Lock-In (demod. number 1) to [amplitude]. Example: lockin.AmplSet 5e-3
lockin.PhaseGet	-	Returns the phase of the Lock-In (demod. number 1). Example: lockin.PhaseGet
lockin.PhaseSet [phase]	[phase]: (°)	Sets the phase of the Lock-In (demod. number 1) to [phase]. Example: lockin.PhaseSet 55.8

Atom tracking

Command	Parameter: (Unit)	Description and Example
atomtrack.ModSet [status]	[status]: (none)	Turns the atom tracking modulation on or off. [status]: 1 – Modulation ON [status]: 0 – Modulation OFF Example: atomtrack.ModSet 1

atomtrack.TrackSet [status]	[status]: (none)	<p>Turns the atom tracking controller on or off.</p> <p>[status]: 1 – Controller ON [status]: 0 – Controller OFF</p> <p>Example: atomtrack.TrackSet 1</p>
------------------------------------	------------------	--

Tip position

Command	Parameter: (Unit)	Description and Example
z.Get	-	<p>Returns the Z position of the tip in meters.</p> <p>Example: z.Get</p>
z.Set [Z]	[Z]: (m)	<p>Sets the Z position of the tip to [Z].</p> <p>Example: z.Set 50e-9</p>
xy.Get	-	<p>Returns the X and Y position of the tip in meters.</p> <p>Example: xy.Get</p>
xy.Set [X] [Y]	[X]: (m) [Y]: (m)	<p>Moves the tip laterally to [X] in x-direction and [Y] in y-direction.</p> <p>Example: xy.Set 5e-9 -30e-9</p>
x.Add [dx]	[dx]: (m)	<p>Moves the tip relative to the current position in x-direction by [dx].</p> <p>Example: x.Add 100e-12</p>
y.Add [dy]	[dy]: (m)	<p>Moves the tip relative to the current position in y-direction by [dy].</p> <p>Example: y.Add -50e-12</p>
z.Add [dz]	[dz]: (m)	<p>Moves the tip relative to the current position in z-direction by [dz].</p> <p>Example: z.Add 10e-12</p>

withdraw	-	Withdraws the tip. Example: withdraw
-----------------	---	---

Drift compensation

Command	Parameter: (Unit)	Description and Example
drift.Get	-	Returns the status and drift compensation values for x-, y- and z-direction. Example: drift.Get
drift.Set [status] [Vx] [Vy] [Vz]	[status]: (none) [Vx]: (m/s) [Vy]: (m/s) [Vz]: (m/s)	Sets the values for the drift compensation in x-, y- and z-direction. [status]: 1 – compensation ON, 0 – compensation OFF [Vx]: x-direction [Vy]: y-direction [Vz]: z-direction Example: drift.Set 1 10e-15 0 30e-15
drift.correctZ [time]	[time]: (s)	Estimates and corrects the drift in z-direction by measuring drift for [time] seconds. Example: drift.correctZ 120

3) Tip Manipulation


Function: Enables manual movement of the tip along the **x**, **y**, and **z** directions.

Usage:

- Set the desired **step size** in the *Step Size* box.
- Move the tip by pressing the directional buttons:
 - **X+ / X-** → Moves the tip in the x-direction.
 - **Y+ / Y-** → Moves the tip in the y-direction.
 - **Z+ / Z-** → Moves the tip in the z-direction.
- Each button press adjusts the tip position by the preset step size.

Important Note:

- Movement in the **z-direction** is only possible when the **feedback loop is switched OFF**.




4) External TCP Interface

Function: Provides an overview of the external TCP interfaces. The TCP interface allows to **send commands to other devices** over a network using the TCP protocol. This enables to mix **local function calls** with **remote device control** inside the same script.

Usage:

- Remote commands are added in the same way as normal commands ([Adding New Functions to the Command System](#)).



5) Adding New Functions to the Command System

Structure of a Registry Entry

The **command system** works by mapping **command names** (the text the user writes in the script) to Python functions or class methods. These can be functions or class methods provided by the **Python Interface Package for Nanonis (nanonis_spm)**, or custom, self-implemented functions for more complex tasks.

This mapping is stored in the **FUNCTION_REGISTRY**, where each entry describes:

- **Which function to call**
- **What arguments it expects**
- **Which arguments are user-provided vs. default**
- **How arguments should be converted (type casting)**

```
FUNCTION_REGISTRY = {
    "command_name": {
        "func": <function_or_method_reference>,
        "args": {
            "arg1_name": {
                "type": <Python_type>,
                "default": <default_value>,
                "user": <True_or_False>
            },
            "arg2_name": {
                "type": <Python_type>,
                "default": <default_value>,
                "user": <True_or_False>
            },
            ...
        }
    }
}
```

Explanation of fields:

command_name:

- The text the user must type in the script to call this function.

func:

- A reference to the function or method to be executed.
- Nanonis commands are made available via the **Python Interface Package for Nanonis**.
- For standalone functions: just use the function name (e.g. complexTask).
- For class methods: use the **unbound method** (e.g. nanonis_spm.Nanonis.Bias_Set).
 - The system will automatically inject the correct instance from INSTANCE_MAP. The INSTANCE_MAP is a lookup table that maps class names to specific pre-created objects, so that the parser knows which instance to use when executing class methods.

args: (dictionary of arguments)

Defines the **arguments** required by the function and their respective **data types**.

Each argument must define:

- **"type":** Python data (numpy) type used to convert the user input (e.g. np.int, np.float32).
- **"default":** Value used if this argument is not provided by the user.
- **"user":**
 - ❖ True → user must provide this argument in the script.
 - ❖ False → argument is filled automatically with the default.

Functions and argument definitions for Nanonis commands follow the **Nanonis TCP protocol**.

Arguments are always read **in the dictionary order**.

Examples

These commands are based on Nanonis functions that are provided via the TCP protocol.

```
"lockin.PhaseSet": {
    "func": nanonis_spm.Nanonis.LockIn_DemodPhasSet,
    "args": {
        "Demodulator number": {"type": np.int32, "default": 1, "user":
False},
        "Phase (deg)": {"type": np.float32, "default": 0.0, "user": True},
    },
},
```

```
"drift.Set": {
    "func": nanonis_spm.Nanonis.Piezo_DriftCompSet,
    "args": {
        "Compensation status": {"type": np.uint32, "default": 1, "user":
True},
        "Vx (m/s)": {"type": np.float32, "default": 0.0, "user": True},
        "Vy (m/s)": {"type": np.float32, "default": 0.0, "user": True},
        "Vz (m/s)": {"type": np.float32, "default": 0.0, "user": True},
        "Saturation limit (%)": {"type": np.float32, "default": 10.0, "user":
False},
    },
}
```

This command is based on a custom, self-implemented function.

```
"bias.Add": {
    "func": ScriptingInterface.addBias,
    "args": {
        "Bias value (V)": {"type": np.float32, "default": 0.1, "user": True},
    },
},
```

Adding a Remote Command

All external devices are defined in the **TCP_INTERFACES** dictionary.

```
TCP_INTERFACES = {  
    "QuPe": {"host": "192.168.1.10", "port": 1337},  
}
```

Multiple TCP interfaces can be managed via **TCP_INTERFACES**. Each remote command specifies which interface to use.

The **TCPClient** class handles communication. It has two methods:

1. **send(command, args)**
 - Sends a command with arguments to the device.
 - Fire-and-forget (does not wait for a response).
2. **query(command, args)**
 - Sends a command with arguments to the device and waits for a response.
 - Expected format from the device: error|response|variable1,variable2
 - "None" for empty values

Remote commands are added to the **FUNCTION_REGISTRY** just like local ones. Functions are implemented as lambda callings.

1. Decide a **command name** (e.g., " QuPe.FreqSet").
2. Create a **lambda function** that calls **TCPClient(interface_name).send(...)**.
3. Define **arguments** in the same style as local commands.

Example:

```
"QuPe.FreqSet": {  
    "func": lambda freq: TCPClient("QuPe").send("setFrequency", [freq]),  
    "args": {  
        "Frequency (Hz)": {"type": np.float32, "default": 1e9, "user": True},  
    },  
}
```

- **freq** is the user-supplied argument.
- Defaults are used if user=False.
- The TCP client sends the command **"setFrequency freq"** (**setFrequency 1e9**) to the interface **QuPe**.